



Pfu DNA polymerase

Cat: PE2010
Conc. 5 units/ μ l

Size: 500 Units
Store at -20°C .

Description

Pfu DNA Polymerase is a thermostable enzyme with a molecular weight of 90 kDa. It catalyzes the polymerization of nucleotides into duplex DNA in the $5' \rightarrow 3'$ direction, resulting in blunt-ended PCR products without 3'-dA overhangs. Pfu DNA Polymerase exhibits $3' \rightarrow 5'$ exonuclease (proofreading) activity that enables the polymerase to correct the mis-incorporation of nucleotide, and lacks $5' \rightarrow 3'$ exonuclease activity. It is suitable for PCR and primer extension reaction that requires high fidelity when the PCR fragment is relatively shorter than 3kb.

The extension rate of Pfu DNA Polymerase is about 600bp/min in standard condition. The appropriate reaction temperature is $65^{\circ}\text{C} \sim 75^{\circ}\text{C}$, the working concentration of dNTPs is 100-300 μ M, the working concentration of Mg^{2+} is 2~3mM, and the optimal pH is 8.1~9.1. The amount of enzyme is 1~1.5unit for 20 μ l PCR reaction, while 2~3units for 50 μ l PCR reaction.

Reaction Buffer (PCR Buffer, 10X):

100mM KCl, 160mM $(\text{NH}_4)_2\text{SO}_4$, 20mM MgSO_4 , 200mM Tris-HCl (pH8.8), 1% Triton X-100, 1mg/ml BSA

Storage Buffer

Pfu DNA Polymerase is supplied in 50 mM Tris-HCl (pH 8.2), 0.1 mM EDTA, 0.1% Tween20, 0.1% NP-40, 1mM DTT, and 50% Glycerol.

Unit Definition

One unit incorporates 10 nanomole of dNTPs into acid insoluble material in 30 minutes at 74°C

Product Contents:

- | | |
|--|---------------------|
| 1. Pfu DNA Polymerase (5 units/ μ l) | 100 μ l (500 U) |
| 2. 10 x PCR reaction buffer | 1.5 ml |

Protocol

- 1- In a 0.2 ml or 0.5 ml thin wall tube, add the following components:

| | 20 μ l reaction | 50 μ l reaction |
|--------------------------|---------------------|---------------------|
| 10 x PCR reaction buffer | 2 μ l | 5 μ l |
| Forward primer | 10~20 pmol | 20~50 pmol |
| Reverse primer | 10~20 pmol | 20~50 pmol |
| Template DNA | 5~50 ng | 10~100 ng |
| dNTPs (10 mM each) | 0.5 μ l | 1 μ l |
| PCR dye (optional) | 2 μ l | 5 μ l |
| Pfu DNA Polymerase | 1-1.5u | 2-3u |
| ddH ₂ O | X μ l | X μ l |
| Total | 20 μ l | 50 μ l |

- 2- Mix gently by pipetting, close the tube and centrifuge for a few seconds.
3- Add mineral oil to each tube (this step is unnecessary when using a thermal cycler with top heating).
4- Perform PCR cycles according to the PCR condition. An example is given as follows:
(Annealing temperature and time need to be optimized for each primer/template combination).

| | | |
|--------------------------|----------|----------------|
| 94 $^{\circ}\text{C}$ | 2.5 min | } 25~35 cycles |
| 94 $^{\circ}\text{C}$ | 45 sec | |
| 50~72 $^{\circ}\text{C}$ | 1 min | |
| 72 $^{\circ}\text{C}$ | 1~3 min | |
| 72 $^{\circ}\text{C}$ | 5~10 min | |

- 5- Run 2~5 μ l PCR products on 1% agarose gel stained with GelRed, GelGreen or Et.Br



Accessory Products

The following products are available separately from Genbiotech.

| Product | Quantity | Catalog no. |
|--------------------------------------|------------|-------------|
| Taq DNA polymerase | 500 U | PE1010 |
| dNTPs in Separate Tube (100 mM each) | 4 x 100 ul | PD2001 |
| dNTPs in Separate Tube (100 mM each) | 4 x 250 ul | PD2002 |
| dNTPs in Separate Tube (100 mM each) | 4 x 400 ul | PD2003 |
| EvaGreen™, 20X in water | 5 x 1 ml | 3100 |
| Loading Dye 6X | 3 x 1 ml. | NM4011 |
| Custom oligonucleotides | | |